

[Home](#) > [Beauty](#)**BEAUTY**

- » Brand directory
- » Beauty Contact List

**HOT LINKS**

- » What's On
- » Mystery Shopping
- » Wedding Shop
- » Gift Card

**EMAIL SIGN UP**

Updates, special offers  
and promotions




JO MALONE


[» Buy Online](#)

BEAUTY

[« back](#)

## Anti-Aging

**STOP™** the clock and reverse signs of ageing with the revolutionary skin renewal device that offers clinical results of skin tightening and reduced fine lines and wrinkles from the comfort of the home.

STOP™ is the perfect choice for those concerned with lined or lax skin wishing to turn back vital signs of ageing without turning to invasive surgery techniques or visiting professional clinics.

STOP™ is an innovative hand held device powered by TriPollar™ Radio Frequency (RF) technology. This technology has already transformed the professional beauty market with anti-ageing treatments that boost skin renewal producing younger, smoother skin. Now the same proven technology is available in a sleek, palm sized device offering clinical results in the home.

Scientifically developed and clinically proven, TriPollar™ technology focuses four low power beams deep into the dermis; much like a laser lens focuses its light. Emitting low power high frequency electrical energy into the dermis, dermal activity is stimulated to tighten collagen fibres and increase new collagen production. Fibroblasts are stimulated increasing collagen replenishment, resulting in dermal thickening. A thicker dermis will make the skin of a woman in her 40's look and feel like she is in her 30's with visible long lasting results - reduction of fine lines and wrinkles and improved skin texture.

From the initial treatment the results of smoother more radiant skin are noticeable. Over a period of 4-6 weeks, repeating the ~15 minute treatment twice per week, long lasting effects are achieved.

Clinical studies showed up to 50% increase in dermal thickening.

The appearances of wrinkles are reduced by 16% and fine lines are reduced by 35% over a 3.5 week STOP™ trial.

### STOP™: technological innovation

Through creating a scientifically-proven product for the home featuring TriPollar™, the third generation of Radio Frequency technology, STOP™ puts the power to control the signs of ageing literally in the hands of the consumer.

### STOP™: how to use

The product, which redefines contours on the face and neck area includes STOP™ Preparation Gel and STOP™ After Treatment.

Simply prepare the face by applying STOP™ Preparation Gel; connect the device to a power supply, select from a setting of low, medium or high. Progress with the treatment by ensuring all four electrodes touch the skin. Glide the device in even, circular motions across the face for approximately 15 minutes and a warm, pleasant, massaging sensation will be felt. Following the treatment apply the STOP™ After Treatment to soothe.

### Time for the clinical way.

Time to STOP™.

#### Site information

- » Privacy
- » Terms and conditions
- » Sitemap

#### Customer services

- » Careers
- » Corporate info
- » Contact us

#### Our shops

- » Store location
- » What's on
- » Contact us

#### About us

- » History
- » FAQs
- » Mystery shopping

#### Useful links

- » Selfridges card
- » Email sign-up
- » Environment